

MANUFACTURED STERILIZED CONTAINERS - ACCEPTANCE TESTS
(Per Lot)

Brand: _____ Lot Number: _____

Volume check: **Bottle type** (circle one): Plain Dechlor

Fluorescence check:

Container type (circle one):

Plain bottle Dechlor bottle

Quanti-Tray 200 Quanti-Tray 2000

Fluorescence present (circle one): Yes / No

Results of fluorescence check are:

acceptable / unacceptable (circle one)

Container	Measured
Number	Volume (ml)
1	

Acceptance Criteria: Volume = 100 mL \pm 2.5 mL = 97.5 mL to 102.5 mL.

Calibration of containers is: *acceptable / unacceptable (circle one)*

Analyst: _____

Analyst: _____

Date: _____

Date: _____

Sterility Check:

Container Type: (circle one): Plain bottle Dechlor bottle Quanti-Tray 200 Quanti-Tray 2000 Whirl-pak bag

One container from (*brand*) _____ lot number _____ was tested for
sterility by aseptically filling it with 25 mL* sterile _____ (non-selective broth).

After 24 hr incubation at 35 \pm 0.5°C, *growth / no growth (circle one)* was observed.

After 48 hr incubation at 35 \pm 0.5°C, *growth / no growth (circle one)* was observed.

Sterility test results for this lot of containers are: *acceptable / unacceptable (circle one)*

Analyst: _____

Date: _____

*Use 100 mL for Quanti-Tray

Sodium Thiosulfate / Chlorine Residual:

Positive Control _____ ppm

Test Container _____ ppm

Acceptance Criteria: Test container - chlorine residual not detected.

Sodium thiosulfate/Chlorine Residual Check is: *acceptable / unacceptable (circle one)*.

Analyst: _____

Date: _____