

PIPETTOR CALIBRATION

Annually

PIPETTOR ID: _____

VOLUME OF PIPETTOR _____

DATE	INITIALS	INDIVIDUAL DELIVERIES g										AVERAGE	
		1	2	3	4	5	6	7	8	9	10		

± 2.5% - acceptable range for water testing

0.1 ml: 0.0975 to 0.1025 g

1 ml: 0.975 to 1.025g

Tolerances based on the weight of water at 19°C to 24°C.

Balance Model # _____ Balance ID# _____ Year: _____
 Weight Set ID# _____ Date of Weight Calibration Check: _____

Date	Low Range		Accept	Mid Range		Accept	High Range		Accept	Initials
	Standard Mass	Balance Reading	Y/N	Standard Mass	Balance Reading	Y/N	Standard Mass	Balance Reading	Y/N	

Comments: _____
